

LaVue: Live viewing images from detectors - supports all 2D detectors used in PETRA III

communication with sardana

Jan Koteński, Christoph Rosemann, André Rothkirch

Deutsches Elektronen-Synchrotron



February 18, 2021

Application Layout and Image Sources

Image source

Image name – Button box

General tools

Preparation:

Background image

Mask image

Transformations


Display intensity:

Scaling

Min./Max. levels

Color gradient

Statistics



2D image – Specialised tools

Application Layout and Image Sources

Image source

Image name – Button box

HIDRA:

Pilatus, Eiger
PCO, Perkin Elemer

HTTP response:

Eiger

Tango Attribute:

Lambda, PCO
Jungfrau, AGIPD
LimaCCDs (e.g. Andor)

Tango Events:

LimaCCDs

Tango File:

Pilatus w/o Hidra


ZMQ Stream:

processed image (P06)

Nexus File:

Nexus Writer (SMWR)

Epics PV, DOOCS, Tine, ASAPO ...



2D image – Specialised tools

Application Layout and Image Sources

Image source

Image name – Button box

In the **expert mode**:


i.e. `lavue -m expert`

General
configuration

Dialog layout

Bookmarks for
image source parameters

Custom color gradients



All stored in `$HOME/.config/DESY/LaVue.conf`

Application Layout and Image Sources

Image source

Image name – Button box

Specialized Tools

Intensity

ROI

LineCut

Angle/Q

MoveMotors

MeshScan

1d-Plot


Projections

Maxima

Parameters

Q+ROI+Proj

Diffraction




Tools added on request

Application Layout and Image Sources

General tools

- Multiple image source
- Range Window and Down-sampling
- Filters
- Memory buffer
- Color/RGB channels
- Image Subtraction
- Image Mask
- Mask high values
- Transformations
- Intensity scale
- Intensity levels
- Histogram
user color gradients
- Statistics



Configuration

General configuration

Image source:

- Number of image sources
- Refresh time in s

Image display:

- Keep original coordinates
- Ranges and ROIs Colors

ROIs:


- Send ROIs values

Sardana:

- Enabled
- Door
- Add ROIs to MG with `nxsadd` macro
- Fetch ROIs order from Sardana Env.

In the **expert mode**:
i.e. `lavue -m expert`

All stored in `$HOME/.config/DESY/LaVue.conf`



Configuration

Layout configuration

Show control bar widget:

- more options


Image Source Selection:

Tool Widget Selection:

- select combobox content
- drag and drop

In the **expert mode**:
i.e. `lavue -m expert`

All stored in `$HOME/.config/DESY/LaVue.conf`



Configuration


Image Sources configuration

Labels for Sources

- star and unstar
- widget with {Label: Value}

In the **expert mode**:
i.e. `lavue -m expert`

All stored in `$HOME/.config/DESY/LaVue.conf`



Configuration

Filter Plugins configuration

Widget with {filter: parameter}

Filters:

- package.module.class
- package.module.function

input: image, imagename,
metadata and image widget
output: processed image


Parameters:

- any string - parameter
of filter class constructor

In the **expert mode**:

i.e. `lavue -m expert -j myinstance`

All stored in `$HOME/.config/DESY/LaVue.conf`



Specialized Tool - Diffractogram

Azimuth integration of defined range

Calibration:
from pyFAI

Range in deg


Units:
 q , 2θ or r

Up to 4
diffractograms

Show and Next
buttons

Results:
sent to tango

LavueController.ToolResults



Lavue Controller - Tango Server

Communicate with LaVue

via Tango interface

- exec: `lavue -a p00/lavuecontroller/1`

Get/Set Detector parameters

- BeamCenterX, BeamCenterY, PixelSizeX, PixelSizeY

- Energy, DetectorDistance

Get/Set ROI bounds

- DetectorROIs, *Values, *Params

Tools: ToolResults

Control LaVue via json LavueState:

import tango


```
lc = tango.DeviceProxy('p09/lavuecontroller/1')
```

```
lc.LavueState = '{"source": "tangoattr", \n    "configuration": "sys/tg_test/1/double_image_ro"}'
```

```
lc.LavueState = '{"start": true}'
```

```
lc.LavueState = '{"tool": "roi"}'
```


```
lc.LavueState = '{"stop": true}'
```



Region Of Interest bounds

Send ROI bounds to a Tango Server by pressing Add button

- exec: `lavue -n p00/lambdaonlineanalysis/1`
- Tango server has a SPECTRUM `RoIs` attribute (python slice convention)



Send ROI bounds to Sardana MacroServer Environment

- variables: `DetectorROIs`, `DetectorROIsOrder`, `DetectorROIsParams`
- ROI alias(es) should be defined (or select 'Single ROIs aliases')

Summary

- Goal: **LaVue** supports all 2D detectors used in PETRA III
- Tuned to DESY needs: we implement **your requests**
- More info:
<https://confluence.desy.de/display/FSEC/LaVue+-+Live+Image+Viewer>
- Sources: <https://github.com/lavue-org/lavue/>
- Debian packages: deb <http://repos.pni-hdri.de/apt/debian>

...

- Goal: **LaVue** supports all 2D detectors used in PETRA III
- Tuned to DESY needs: we implement **your requests**
- More info:
<https://confluence.desy.de/display/FSEC/LaVue+-+Live+Image+Viewer>
- Sources: <https://github.com/lavue-org/lavue/>
- Debian packages: deb <http://repos.pni-hdri.de/apt/debian>

Thank You !